



Recombinant Human Interferon regulatory factor 2-binding protein-like (IRF2BPL)

Product Code	CSB-MP887953HU
Abbreviation	IRF2BPL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H1B7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSAAQVSSSR RQSCYLCDLP RMPWAMIWDF SEPVCRGCVN YEGADRIEFV IETARQLKRA HGCFQDGRSP GPPPPVGVKT VALSAKEAAA AAAAAAAAAA AAQQQQQQQQ QQQQQQQQQQ QQQQQQQLNH VDGSSKPAVL AAPSGLERYG LSAAAAAAAA AAAAVEQRSR FEYPPPPVSL GSSSHTARLP NGLGGPNGFP KPTPEEGPPE LNRQSPNSSF AAASVASRRG THGGLVTGLP NPGGGGGPQL TVPPNLLPQT LLNGPASAASV LPPPPPHALG SRGPPTPAPP GAPGGPACLG GTPGVSAATSS SASSTSSSV AEVGVGAGGK RPGSVSSTDQ ERELKEKQRN AEALAESES LRNRREEWAS KPKMVRDTLL TLAGCTPYEV RFKKDHSLLG RVFAFDAVSK PGMDYELKLF IEYPTGSGNV YSSASGVAKQ MYQDCMKDFG RGLSSGFKYL EYEKKHGSGD WRLLDLLPE AVRFFKEGVP GADMLPQPYL DASCPLPTA LVSLSRAPSA PPGTGALPPA APSGRGAAAS LRKRKASPEP PDSAEGALKL GEEQQRQQWM ANQSEALKLT MSAGGFAAPG HAAGGPPPPP PPLGPHSNRT TPESAPQNG PSPMAALMSV ADTLGTAHSP KDGSSVHSTT ASARRNSSSP VSPASVPGQR RLASRNGDLN LQVAPPPPSA HPGMDQVHPQ NIPDSPMANS GPLCCTICHE RLEDTHFVQC PSVPSHKFCF PCSRESIKAQ GATGEVYCPS GEKCPLVGSN VPWAFMQGEI ATILAGDVKV KKERDP
Source	Mammalian cell
Target Names	IRF2BPL
Protein Names	Recommended name: Interferon regulatory factor 2-binding protein-like Alternative name(s): Enhanced at puberty protein 1
Expression Region	1-796
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.