



Recombinant Human Neurogenin-2 (NEUROG2)

Product Code	CSB-BP887960HU
Abbreviation	NEUROG2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H2A3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFVKSETLEL KEEEDVLVLL GSASPALAAL TPLSSSADEE EEEEPGASGG ARRQRGAEEAG QGARGGVAAG AEGCRPARLL GLVHDCKRRP SRRAVSRGA KTAETVQRIK KTRRLKANNR ERNRMHNLNA ALDALREVLV TFPEDAKLTK IETLRFAHNY IWALTETLRL ADHCGGGGGG LPGALFSEAV LLSPGGASAA LSSSGDSPSP ASTWSCTNSP APSSSVSSNS TSPYSCTLSP ASPAGSDMDY WQPPPPDKHR YAPHLPIARD CI
Source	Baculovirus
Target Names	NEUROG2
Protein Names	Recommended name: Neurogenin-2 Short name= NGN-2 Alternative name(s): Class A basic helix-loop-helix protein 8 Short name= bHLHa8 Protein atonal homolog 4
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.