



Recombinant Human Uncharacterized protein C9orf37 (C9orf37)

Product Code	CSB-MP887963HU
Abbreviation	C9orf37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H2J1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEFHYVAQAD LELLTSSNPP ASASQSTGIT GGSHRARP GP VHFIDKVTDK PSHSHPFALK ENWNLNPEPS SPPSPLFLEA PSRQASQHHG ASPGAGTSAG CPFEKCCSTE PCLSGLGDVG RGEAASLRAR PGSGASRGQG PGSRVSCRRD LGKPLHAPAG FSAGEVHTTP LGNLGA
Source	Mammalian cell
Target Names	ARRDC1-AS1
Protein Names	Recommended name: Uncharacterized protein C9orf37
Expression Region	1-176
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	While this gene appears to be transcribed and translated to make a protein, no function has been assigned to the protein.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.