



Recombinant Human PR domain zinc finger protein 13 (PRDM13)

Product Code	CSB-EP887981HU-B
Abbreviation	PRDM13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H4Q3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MHGAARAPAT SVSADCCIPA GLRLGPVPGT FKLKGYLSDR REPGPKKKVR MVRGELVDES GGSPLEWIGL IRAARNSQEQ TLEAIADLPG GQIFYRALRD VQPGEELTVW YSNSLAQWFD IPTTATPTH D EKGEERYICW YCWRTFRYPN SLKAHLRFHC VFSGGGGGAF LHHEHAARQG AVPAADGLGL SPKPPAPDFA APSQAGTLRP HPLGPPPVA CGAREGIKRE ASSAPSATSP TPGKWWGQPKK GKEQLDRALD MSGAARGQGH FLGIVGGSSA GVGLAFYYPG VRSAFKPAGL ARAAAAHGD PYREESSKQ GAGLALGRLL GGGRACGRPG SGENSAAGGA GHHHHHHAHH HHPKCLLAG DPPPPPPGL PCSGALRGFP LLSVPPEEAS AFKHVERAPP AAAALPGARY AQLPPAPGLP LERCALPLD PGGLKAYPGG ECSHLPVAVMP AFTVYNGELL YGSPATTAYY PLKLHFGGLL KYPESISYFS GPAAAALSPA ELGSLASIDR EIAMHNQQLS EMAAGKGRGR LDSGTLPPAV AAAGGTGGGG SGGSGAGKPK TGHLCLYCGK LYSRKYGLKI HMRTHTYGYP LKCKVCLRPF GDPSNLNKHI RLHAEGNTPY RCEFCGKVLV RRRDLERHVK SRHPGQSLLA KAGDGPGAEP GYPPEPGDPK SDDSDVDVCF TDDQSDPEVG GGERDL
Source	E.coli
Target Names	PRDM13
Protein Names	Recommended name: PR domain zinc finger protein 13 Alternative name(s): PR domain-containing protein 13
Expression Region	1-707
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.