



Recombinant Human Ribonuclease P protein subunit p21 (RPP21)

Product Code	CSB-MP887986HU
Abbreviation	RPP21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H633
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AGPVKDREA FQRLNFLYQA AHCVLAQDPE NQALARFYCY TERTIAKRLV LRRDPSVKRT LCRGCSSLLV PGLTCTQRQR RCRGQRWTVQ TCLTCQRSQR FLNDPGHLLW GDRPEAQLGS QADSKPLQPL PNTAHSISDR LPEEKMQTQG SSNQ
Source	Mammalian cell
Target Names	RPP21
Protein Names	Recommended name: Ribonuclease P protein subunit p21 Short name= RNaseP protein p21 EC= 3.1.26.5 Alternative name(s): Ribonuclease P/MRP 21 kDa subunit Ribonucleoprotein V
Expression Region	2-154
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.