



Recombinant Human DENN domain-containing protein 2D (DENND2D)

Product Code	CSB-MP887987HU
Abbreviation	DENND2D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6A0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEGQVVGRVF RLFQRLLQL RAGPPQDNSG EALKEPERAQ EHSLPNFAGG QHFFEYLLVV SLKKKRSEDD YEPITYQFP KRENLLRGQQ EEEERLLKAI PLFCFPDGNE WASLTEYPRE TFSFVLTNVD GSRKIGYCRR LLPAGPGPRL PKVYCIISCI GCFGLFSKIL DEVEKRHQIS MAVIYPFMQG LREAAFPAPG KTVTLKSFIP DSGTEFISLT RPLDSHLEHV DFSSLLHCLS FEQILQIFAS AVLERKIIFL AEGLSTLSQC IHAAAALLYP FSWAHTYIPV VPESLLATVC CPTPFMVG VQ MRFQQEVMDS PMEEVLLVNL CEGTFLMSVG DEKDILPPKL QDDILD SLGQ GINELKTAEQ INEHVSGPFV QFFVKIVGHY ASYIKREANG QG HFQERSFC KALTSKTNRR FVKKFVK TQL FSLFIQEA EK SKNPPAGYFQ QKILEYEEQK KQKKPREKTV K
Source	Mammalian cell
Target Names	DENND2D
Protein Names	Recommended name: DENN domain-containing protein 2D
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.