



Recombinant Human Transcription factor SOX-17 (SOX17)

Product Code	CSB-EP887989HU-B
Abbreviation	SOX17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6I2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSSPDAGYAS DDQSQTQSAL PAVMAGLGPC PWAESLSPIG DMKVKGEAPA NSGAPAGAAG RAKGESRIRR PMNAFMVWAK DERKRLAQQN PDLHNAELSK MLGKSWKALT LAEKRPFVEE AERLRVQHMQ DHPNYKYRPR RRKQVKRLKR VEGGFLHGLA EPQAAALGPE GGRVAMDGLG LQFPEQGFP A GPPLLPPHMG GHYRDCQSLG APPLDGYPLP TPDTSPLDGV DPDPAFFAAP MPGD CPAAGT YSYAQVSDYA GPPEPPAGPM HPRLGPEPAG PSIPGLLAPP SALHVYYGAM GSPGAGGGRG FQMQPQHQQH HQQHQQHPPGP GQPSPPPEAL PCR DGT DPSQ PAELLGEVDR TEFEQYLHFV CKPEMGLPYQ GHDSGVNLPD SHGAISSVVS DASSAVYYCN YPDV
Source	E.coli
Target Names	SOX17
Protein Names	Recommended name: Transcription factor SOX-17
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.