



Recombinant Human Lysine-specific demethylase NO66 (NO66)

Product Code	CSB-EP887992HU-B
Abbreviation	NO66
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6W3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDGLQASAGP LRRGRP KRRR KPQPHSGSVL ALPLRSRKIR KQLRSVVSRM AALRTQTLPS ENSEESRVES TADDLGDALP GGAAVA AAVPD AARREPYGHL GPAELLEASP AARSLQTPSA RLVPASAPPA RLVEVPAAPV RVVETSALLC TAQHLLAAVQS SGAPATASGP QVDNTGGGEP WDSPLRRVLA ELNRIPSSRR RAARLFEWLI APMPPDHFYR RLWEREAVLV RRQDHTYYQG LFSTADLDSM LRNEEVQFGQ HLDAARYING RRET LNPPGR ALPAAAWSLY QAGCSLRLLC PQAFSTTVWQ FLAVLQEQFG SMAGSNVYLT PPNSQGFAPH YDDIEAFVLQ LEGRKLWRVY RPRVPTTEELA LTSSPNFSQD DLGEPVLQTV LEPGDLLYFP RGGFIHQAEQ DGVHSLHLTL STYQRNTWGD FLEAILPLAV QAAMEENVEF RRGLPRDFMD YMGAQHSDSK DPRRTAFMEK VRVLVARLGH FAPVDAVADQ RAKDFIHDSL PPVLTDRERA LSVYGLPIRW EAGEPVNVGA QLTTETEVMH LQDGIARLVG EGGHFLYYT VENS RVYHLE EPKCLEIYPQ QADAMELLLG SYPEFVRVGD LPCDSVEDQL SLATTLYDKG LLLTKMPLAL N
Source	E.coli
Target Names	RIOX1
Protein Names	Recommended name: Lysine-specific demethylase NO66 EC= 1.14.11.27 Alternative name(s): Nucleolar protein 66 Short name= hsNO66
Expression Region	1-641
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.