



Recombinant Human UPF0672 protein CXorf36 (CXorf36)

Product Code	CSB-EP887998HU-B
Abbreviation	CXorf36
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H7Y0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	LPASSLSSL VPQVRTSYNF GRTFLGLDKC NACIGTSICK KFFKEEIRSD NWLASHLGLP PDSLLSY PAN YSDDSKIWRP VEIFRLVSKY QNEISDRRIC ASASAPKTCS IERVLRKTER FQKWLQAKRL TPDLVQGLAS PLLRCPSQRL LDRVVRRYAE VADAGSIFMD HFTDRDKLRL LYTLAVNSHP ILLQIFPGAE GWPLPKYLGS CGRFLVSTST RPLQEFYDAP PDQAADLAYQ LLGVLESLRS NDLNYFFYFT HIDAGMFGVF NNGHLFIRDA SAVGVIDKQE GSQEANRAGE NKDIFSCLVS GCQAQLPSCE SISEKQSLVL VCQKLLPRL QGRFPSPVQD DIDSILVQCG DSIRPDPEVL GAASQLKDIL RPLRTCDSRF AYRYPDCKYN DKF
Source	E.coli
Target Names	CXorf36
Protein Names	Recommended name: UPF0672 protein CXorf36
Expression Region	32-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.