



# Recombinant Human Bromodomain-containing protein 9 (BRD9)

<b>Product Code</b>	CSB-BP888000HU
<b>Abbreviation</b>	BRD9
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H8M2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGKKHKHKA EWRSSYEDYA DKPLEKPLKL VLKVGGSSEVT ELSGSGHDSS YYDDRSDDHER ERHKEKKKKK KKKSEKEKHL DDEERRKRKE EKRRKREREH CDTEGEADDF DPGKKVEVEP PPDRPVRACR TQPAENESTP IQQLEHFLR QLQRKDPHGF FAFPVTDAIA PGYSMIKHP MDFGTMKDKI VANEYKSVTE FKADFKLMCD NAMTYNRPDT VYYKLAKKIL HAGFKMMSKQ AALLGNEDTA VEEPVPEVVP VQVETAKKSK KPSREVISCM FEPEGNACSL TDSTAEHV ALVEHADEA RDRINRFLPG GKMGYLKRNG DGSLLYSVNVN TAEPDADEEE THPVDLSSLS SKLLPGFTTL GFKDERRNKV TFLSSATTAL SMQNNSVFGD LKSDMELLY SAYGDETVQ CALSLQEFVK DAGSYSKKVV DDLDDQITGG DHSRTLFLQK QRRNVPMKPP DEAKVGDVTLG DSSSSVLEFM SMKSYPDVSV DISMLSSLGK VKKELDPDDS HLNLDDETTKL LQDLHEAQAE RGGSRPSSNL SSLSNASERD QHHLGSPSRL SVGEQPDVTH DPYEFLQSPE PAASAKT
<b>Source</b>	Baculovirus
<b>Target Names</b>	BRD9
<b>Protein Names</b>	Recommended name: Bromodomain-containing protein 9 Alternative name(s): Rhabdomyosarcoma antigen MU-RMS-40.8
<b>Expression Region</b>	1-597
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.