



Recombinant Human Bromodomain-containing protein 9 (BRD9)

Product Code	CSB-EP888000HU
Abbreviation	BRD9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H8M2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGKKHKKHKA EWRSSYEDYA DKPLEKPLKL VLKVGGSSEVT ELSGSGHDSS YYDDRSDDHER ERHKEKKKKK KKKSEKEKHL DDEERRKRKE EKRRKREREH CDTEGEADDF DPGKKVEVEP PPDRPVRACR TQPAENESTP IQQLLEHFLR QLQRKDPHGF FAFPVTDAIA PGYSMIKHP MDFGTMKDKI VANEYKSVTE FKADFKLMCD NAMTYNRPDT VYYKLAKKIL HAGFKMMSKQ AALLGNEDTA VEEPVPEVVP VQVETAKKSK KPSREVISCM FEPEGNACSL TDSTAEHV ALVEHADEA RDRINRFLPG GKMGYLKRNG DGSLLYSVNVN TAEPDADEEE THPVDLSSLS SKLLPGFTTL GFKDERRNKV TFLSSATTAL SMQNNVDFGD LKSDMELLY SAYGDETVGQ CALSLQEFVK DAGSYSKKVV DDLDDQITGG DHSRTLFLQK QRRNVPMKPP DEAKVGDVTLG DSSSSVLEFM SMKSYPDVSV DISMLSSLGK VKKELDPDDS HLNLDDETTKL LQDLHEAQAE RGGSRPSSNL SSLSNASERD QHHLGSPSRL SVGEQPDVTH DPYEFLQSPE PAASAKT
Source	E.coli
Target Names	BRD9
Protein Names	Recommended name: Bromodomain-containing protein 9 Alternative name(s): Rhabdomyosarcoma antigen MU-RMS-40.8
Expression Region	1-597
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.