



Recombinant Human WD repeat-containing protein 76 (WDR76)

Product Code	CSB-MP888003HU
Abbreviation	WDR76
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H967
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSRSGAAAEK ADSRQRPQMK VNEYKENQNI AYWVSLRPAQT TVLIKTAKVY LAPFSLSNYQ LDQLMCPKSL SEKNSNNEVA CKKTKIKKTC RRIIPPKMKN TSSKAESTLQ NSSSAVHTES NKLQPKRTAD AMNLSVDVES SQDGDSEDET TPSLDFSGLS PYERKRLKNI SENADFFASL QLSESAARLR EMIEKRQPPK SKRKKPKREN GIGCRRSMRL LKVDPSGVSL PAAPTPTLV ADETPLLPPG PLEMTSENQE DNNERFKGFL HTWAGMSKPS SKNTEKGLSS IKSYPKANLNG MVIS EDTVYK VTTGPIFSMA LHPSETRTL AVGAKFGQVG LCDLTQQPKE DGVYVFHPS QPVSCLYFSP ANPAHILSLS YDGTLCRGDF SRAIFEEVYR NERSSFSSFD FLAEDASTLI VGHWDGNMSL VDRRTPGTSY EKLTS SSMGK IRTVHVHPVH RQYFITAGLR DTHIYDARRL NSRRSQPLIS LTEHTKSIAS AYFSPLTGNR VVTTADCNL RIFDSSCISS KIPLLTIRH NTFTGRWLTR FQAMWDPKQE DCVIVGSMAS PRRVEIFHET GKRVSFSGGE YLVSVC SINA MHPTRYILAG GNSSGKIHFV MNEKSC
Source	Mammalian cell
Target Names	WDR76
Protein Names	Recommended name: WD repeat-containing protein 76
Expression Region	1-626
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.