



Recombinant Human WD repeat-containing protein 76 (WDR76)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-YP888003HU |
| Abbreviation | WDR76 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9H967 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MSRSGAAAEK ADSRQRPQMK VNEYKENQNI AYWVSLRPAQT TVLIKTAKVY LAPFSLSNYQ LDQLMCPKSL SEKNSNNEVA CKKTKIKKTC RRIIPPKMKN TSSKAESTLQ NSSSAVHTES NKLQPKRTAD AMNLSVDVES SQDGDSEDET TPSLDFSGLS PYERKRLKNI SENADFFASL QLSESAARLR EMIEKRQPPK SKRKKPKREN GIGCRRSMRL LKVDPSGVSL PAAPTPTLV ADETPLLPPG PLEMTSENQE DNNERFKGFL HTWAGMSKPS SKNTEKGLSS IKSYPKANLNG MVIS EDTVYK VTTGPIFSMA LHPSETRTL AVGAKFGQVG LCDLTQQPKE DGVYVFHPHS QPVSCLYFSP ANPAHILSLS YDGTLCRGDF SRAIFEEVYR NERSSFSSFD FLAEDASTLI VGHWDGNMSL VDRRTPGTSY EKLTSSSMGK IRTVHVHPVH RQYFITAGLR DTHIYDARRL NSRRSQPLIS LTEHTKSIAS AYFSPLTGNR VVTTADCANL RIFDSSCISS KIPLLTIRH NTFTGRWLTR FQAMWDPKQE DCVIVGSMAS PRRVEIFHET GKRVSFSGGE YLVSVCISINA MHPTRYILAG GNSSGKIHFV MNEKSC |
| Source | Yeast |
| Target Names | WDR76 |
| Protein Names | Recommended name: WD repeat-containing protein 76 |
| Expression Region | 1-626 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.