



Recombinant Human Histone deacetylase complex subunit SAP30L (SAP30L)

Product Code	CSB-BP888011HU
Abbreviation	SAP30L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HAJ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNGFSTEEDS REGPPAAPAA AAPGYGQSCC LIEDGERCVR PAGNASFSKR VQKSISQKKL KLDIDKSVRH LYICDFHKNF IQSVRNKRKR KTSDDGGDSP EHDTDIPEVD LFQLQVNTLR RYKRHYKLQT RPGFNKAQLA ETVSRHFRNI PVNEKETLAY FIYMVKS NKS RLDQKSEGGK QLE
Source	Baculovirus
Target Names	SAP30L
Protein Names	Recommended name: Histone deacetylase complex subunit SAP30L Alternative name(s): HCV non-structural protein 4A-transactivated protein 2 Sin3 corepressor complex subunit SAP30L Sin3-associated protein p30-like
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.