



Recombinant Human UPF0160 protein MYG1, mitochondrial (C12orf10)

Product Code	CSB-EP888015HU-B
Abbreviation	C12orf10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HB07
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GTH NGTFHCDEAL ACALLRLLPE YRDAEIVRTR DPEKCLASCDI VVDVGGEYDP RRHRYDHHQR SFTETMSSL S PGKPWQTKLS SAGLIYLFHG HKLLAQLLGT SEEDSMVGTLYDKMYENFVE EVDVDNGIS QWAEGEPRYA LTTTLSARVA RLNPTWNHPD QDTEAGFKRA MDLVQEEFLQ RLDFYQHSWL PARALVEEAL AQRQVDPSPG EIVELAKGAC PWKEHLYHLE SGLSPPVAIF FVIYTDQAGQ WRIQCVPKEP HSFQSRLPLP EPWRGLRDEA LDQVSGIPGC IFVHASGFTG GHHTREGALS MARATLAQRS YLPQIS
Source	E.coli
Target Names	C12orf10
Protein Names	Recommended name: UPF0160 protein MYG1, mitochondrial
Expression Region	48-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.