



Recombinant Schizosaccharomyces pombe

Uncharacterized protein PB1A10.02

(SPAPB1A10.02)

Product Code	CSB-MP888040SXV
Abbreviation	SPAPB1A10.02
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HDY7
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MCNRQPKAVD LPPNYSCPHL LSFQSVEFNT NPTACDDVFC KRIESEKKYN DFLESLFKKY GRDTSDIAD E VDLATGEIIV NNGHLEALKT KDDIWDPTFN NLEISASNGY EKKLDSSIGN PGEKAVSPVH IEDFQSPQIY KFKNLSLRDE MVSDCVFADE VPLASLFVEN VCNETIPSQS CVRLKINDKT RKVDASALEK KSCLLPNSSG TLTDQRGLDT IKHKSIEQNE ILHVISDTLS SPRRRNPLLS SPKTPLRRSF SKSKVRNSNS TKRRNFISLI SMISPRPNLS THHFNLGFQP LSQQTFSFGS STQNPSSST CKKAFCFQCI SESKKC
Source	Mammalian cell
Target Names	SPAPB1A10.02
Protein Names	Recommended name: Uncharacterized protein PB1A10.02
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.