



Recombinant Halobacterium salinarum Peptide chain release factor subunit 1 (prf1)

Product Code	CSB-MP888112HTM
Abbreviation	prf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HNF0
Product Type	Recombinant Protein
Immunogen Species	Halobacterium salinarum (strain ATCC 700922 / JCM 11081 / NRC-1) (Halobacterium halobium)
Purity	>85% (SDS-PAGE)
Sequence	MSEQDEVPSE DRRKYEFRKV IEELKDYEGS GTQLVTIYIP PDKQISDVVA HVTQEHESEAS NIKSKQTRTN VQDALTSIKD RLRYYDTFPP DNGMVVFSGA VDSGGGRTDM VTEVLESPPQ PIESFRYHCD SAFLTEPLAE MLGDKGLYGL IVLDRRESNV GWLKGKRVQP VKSAESLVPK KQRKGGQSAQ RFARLRLEAI DNFYQEVAGM ADDLFVPKRH EIDGILVGGP SPTKDEFDLDG DYLLHHELQDK VLGKFDVSYT DESGLSDLVD AGQAALAEAD LMDDKSDMEE FFEELNGGKL ATYGFEQTRR NLIMGSVDR LVS EDLRE DV VIYEC PNDHE EYETIDRRNT SPEHTCSDCG EEATEVDRED AIDHLMSIAD QRG TETHFIS TDFEKGEQLL TAFGGYAGIL RYSTGV
Source	Mammalian cell
Target Names	prf1
Protein Names	Recommended name: Peptide chain release factor subunit 1 Alternative name(s): Translation termination factor aRF1
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.