



Recombinant Mouse Pyrin (Mefv), partial

Product Code	CSB-EP888319MO-B
Abbreviation	Mefv
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9JJ26
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	QSLPQCPRHMKQVLLLLFCEDHREPICLICRLSLEHQGHRVVRPIEEALEYKEQI REQLERLREMRGYVEEHRLQGDKKTDDFLKQTEIQKQKISCPLEKLYQLLEKQ EQLFVTWLQELSQTISKVRETYTRVSLLEMIIELEAKQDQPEWDLMQDIGIT LHRAKMMASSELLDTPPGVKEKLHLLYQKSKSVEKNMQCFSEMLSSSEMAFSA SDVAKWEGRQPSATQVQGLVPTVHLKCDGAHTQDCDVVFYPEREAGGSEPK DYLHPQPAQDTPELHEIHSRNNKRKFKSFLKWKPSFSRTDWRLRTCCYRDL QAAAHPNLIFSMI
Source	E.coli
Target Names	Mefv
Protein Names	Recommended name: Pyrin Alternative name(s): Marenosttrin
Expression Region	439-767
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Target Details	This gene encodes a protein, also known as pyrin or marenosttrin, that is an important modulator of innate immunity. Mutations in this gene are associated with Mediterranean fever, a hereditary periodic fever syndrome.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.