



Recombinant Arabidopsis thaliana Probable monodehydroascorbate reductase, cytoplasmic isoform 3 (At3g52880)

Product Code	CSB-EP888714DOA
Abbreviation	At3g52880
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LFA3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAEKSFKYII LGGGVSAGYA AKEFANQGVQ PGELAVISKE AVAPYERPAL SKGYLPEGA ARLPGFHCCV GSGGEKLLPE SYKQKGIELI LSTEIVKADL SAKSLVSATG DVFKYQTLII ATGSTVLRLT DFGVKGADSK NILYLREIDD ADKLVEAIKA KKGKAVVVG GGYIGLELSA VLRINNLDVT MVFPEPWCMF RLFTADIAAF YETYYTNKGV KIIKGTVASG FTAQPNGEVK EVQLKDGRTL EADIVIVGVG AKPLTSLFKG QVEEDKGGIK TDAFFKTSVP DVYAVGDVAT FPLKMYGDVR RVEHVDHSRK SAEQAVKAIK AAEGGA AVEE YDYLPPFFYSR SFDLSWQFYG DNVGDSVLFV DSNPSNPKPR FGAYWVQGGK VVGAFMEGGS GDENKALAKV AKARPSAESL DELVKQGISF AAKI
Source	E.coli
Target Names	MDAR1
Protein Names	Recommended name: Probable monodehydroascorbate reductase, cytoplasmic isoform 3 Short name= MDAR 3 EC= 1.6.5.4
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.