



Recombinant Arabidopsis thaliana E2F transcription factor-like E2FD (E2FD)

Product Code	CSB-YP888719DOA
Abbreviation	E2FD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LFQ9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDSLALAPQV YSRKDKSLGV LVANFLTLYN RPDVDLFGLD DAAAKLGVER RRIYDVVNIL ESIGLVARSG KNQYSWKGFG AVPRALSELK EEGMKEKFAI VPFVAKSEMV VYEKEGEESF MLSPDDQEFS PSPRPDNRKE RTLWLLAQNF VKLFLCSDDD LVTFDSATKA LLNESQDMNM RKKVRRLYDI ANVFSSMKLI EKTHVPETKK PAYRWLGSKT IFENRFIDGS ASLCDRNVPK KRAFGTELTN VNAKRKSGC SKEDSKRNGN QNTSIVIKQE QCDDVKPDVK NFASGSSTPA GTSESNDMGN NIRPRGRLGV IEALSTLYQP SYCNPPELLGL FAHYNETFRS YQEEFGREK
Source	Yeast
Target Names	E2FD
Protein Names	Recommended name: E2F transcription factor-like E2FD Alternative name(s): DP-E2F-like protein 2 E2F-like repressor E2L1
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.