



Recombinant Arabidopsis thaliana Putative F-box protein At5g15660 (At5g15660)

Product Code	CSB-EP888721DOA
Abbreviation	At5g15660
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LFW0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRRRSKKIKT ENNSNPETSE ERNKFDEIPH DLVIEILERL PLKSVARFLT VSKLWATTIR SPDFRKSRYG GSSSEPRTLI VSDLNFKEPN PKLHFFRPSI SSPSFLSSLT CPFTYPRHEE YYYHHVNGLI SVGYGTDQIV INPTTGKFIT LPRPKTRRKL VISFFGYDSV SDQYKVL CMT ERLRGHP EEA SSQH QVYTLG AKQKSWKMIN CSIPHRPWSW NAVCINGVVY YIAKTGEGMF RRCLMRFDLK SDNLDLCTIL PEEIQTSLHD YFLINYKGV AIPNQPNFYT YDVWVMNQEG GKIEWLKNIT FTIKPRKGFV RYLFVTGTTH TGEFILAPTS YTDEFYVFHY NPDMSFRKI RVQAPGVKFS FAQKASVVFS DHSLLRLDNL HIRGSTHTAT GEFILAPRFY SDDLNVIHFN PDTNSFRSTK VEVYEDYE
Source	E.coli
Target Names	At5g15660
Protein Names	Recommended name: Putative F-box protein At5g15660
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.