



# Recombinant Arabidopsis thaliana Putative F-box protein At5g15660 (At5g15660)

<b>Product Code</b>	CSB-MP888721DOA
<b>Abbreviation</b>	At5g15660
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LFW0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRRRSKKIKT ENNSNPETSE ERNKFDEIPH DLVIEILERL PLKSVARFLT VSKLWATTIR SPDFRKSRYG GSSSEPRTL I VSDLNFKEPN PKLHFFRPSI SSPSFLSSLT CPFTYPRHEE YYYHHVNGLI SVGYGTDQIV INPTTGKFIT LPRPKTRRKL VISFFGYDSV SDQYKVL CMT ERLRGHP EEA SSQH QVYTLG AKQKSWKMIN CSIPHRPWSW NAVCINGVVY YIAKTGEGMF RRCLMRFDLK SDNLDLCTIL PEEIQTSLHD YFLINYKGV AIPNQPNFYT YDVWVMNQEG GKIEWLKNIT FTIKPRKGFV RYLFVTGTTH TGEFILAPTS YTDEFYVFHY NPDMSFRKI RVQAPGVKFS FAQKASV VFS DHSLLRLDNL HIRGSTHTAT GEFILAPRFY SDDLNV IHFN PDTNSFRSTK VEVYEDYE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At5g15660
<b>Protein Names</b>	Recommended name: Putative F-box protein At5g15660
<b>Expression Region</b>	1-438
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.