



Recombinant Arabidopsis thaliana F-box/FBD/LRR-repeat protein At3g26920 (At3g26920)

Product Code	CSB-MP888736DOA
Abbreviation	At3g26920
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LJF9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTMKRCKSCG QLSNEDRIS QLPEALLLQI LLLPTKEVV AVSVLAKRWR FLWKMVPSLE FFYYFTNDLE RFSYNSKCL FSHQAPFLQS LHLNMNFGCD PRIMDFEILI GIAFGRQLRK LVLKVYSGDW FKFPSTLYNS ETLETLELYH CILIDVPPV CLKSLRTLNL HEVEFVNDES VVNLLAGCIS LENLVIHQTT DLNVKTFTIA VPSLQRLTVI VEYEEFSVF VVNTPSLKYL KIEGIIVDDR TCIIENTPEL VEASIIDVSF KVFESILGSL ASVQRLSLKV SLVEIFSLPP ISNTFYHLTY LELSTYKPKW WNLLTLMLDT SPNLQVLKIF DFMTSQEQRP WEKWNEPKNV PECLLLHLET FWWTCEYEGKL ENEIELAKYI LRNARRLKKA TFSIIEINPD KRVEVMGELK SVVRASNQCQ LVFI
Source	Mammalian cell
Target Names	At3g26920
Protein Names	Recommended name: F-box/FBD/LRR-repeat protein At3g26920
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.