



Recombinant Arabidopsis thaliana Probable protein phosphatase 2C 13 (At1g48040)

Product Code	CSB-EP888763DOA-B
Abbreviation	At1g48040
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LNF4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MILSQTMVAE AEIRVLDVKC HISAPKDQKN FQIDEVRVSE SVRAEISGSA ETPRFGSGMS CVTTTIGESA SDFIPTIRSG SFADIRSRET MEDEHICIDD LSAHLGSYNF SVPSAFYGVF DGHGGPEAAI FMKENLTRLF FQDAVFPEMP SIVDAFFLEE LENSHRKAFA LADLAMADET IVSGSCGTTA LTALIIGRHL LVANAGDCRA VLCRRGVAVD MSFDHRSTYE PERRRIEDLG GYFEDGYLNG VLAVTRAIGD WELKNPFTDS SSPLISDPEI GQIILTEDDE FLILACDGIW DVLSSQNAVS NVRQGLRRHG DPRQCAMELG KEAARLQSSD NMTVIVICFS SVPSSPKQPQ RRRLRFCVSD EARARLQAML AGE
Source	E.coli
Target Names	At1g48040
Protein Names	Recommended name: Probable protein phosphatase 2C 13 Short name= AtPP2C13 EC= 3.1.3.16
Expression Region	1-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.