



Recombinant Arabidopsis thaliana Allene oxide cyclase 2, chloroplastic (AOC2)

Product Code	CSB-BP888793DOA
Abbreviation	AOC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LS02
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	PSK VQELSVYEIN ELDRHSPKIL KNAFSLMFGL GDLVPFTNKL YTGDLKRRVG ITAGLCVVIE HVPEKKGERF EATYSFYFGD YGHLSVQGPY LTYEDSFLAI TGGAGIFE GA YGQVKLQQLV YPTKLFYTFY LKGLANDLPL ELTGTPVPPS KDIEPAPEAK ALEPSGVISN YTN
Source	Baculovirus
Target Names	AOC2
Protein Names	Recommended name: Allene oxide cyclase 2, chloroplastic EC= 5.3.99.6
Expression Region	78-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.