



Recombinant *Oryza sativa* subsp. japonica Proteasome subunit beta type-2 (PBD1)

| | |
|--------------------------|---|
| Product Code | CSB-EP888804OFG |
| Abbreviation | PBD1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9LST6 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. japonica (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MECVIGVVGR DFAVVAADTS AVQSILVHKT DEDKVMVLDS HKLMGASGEP GDRVQFTEFI QKNLHLYQFR NNIPLSTAAT ANFTRGELAT ALRKNPYYVN VLLAGYDSDV GASLYYIDYI ATFHKIEKGA FGYGSYFCLS LMDKLYRPDM SVEEAVDLVD KCIKEIRLRL VVAPQNFIIK IVDKEGAREY ARRAYTDSPP EAATSEAATV AA |
| Source | E.coli |
| Target Names | PBD1 |
| Protein Names | Recommended name: Proteasome subunit beta type-2 EC= 3.4.25.1 Alternative name(s): 20S proteasome alpha subunit D 20S proteasome subunit beta-4 |
| Expression Region | 1-212 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |