



Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 39 (At3g19360)

Product Code	CSB-MP888812DOA
Abbreviation	At3g19360
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LT81
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDSSYSDSRP MFVQSPYGGW NQTQMIDSMA NPMNNEQGDL HSLSESQSQS QPSQQLQPAL KRPRLVDDNL FNPASSFPQP SSSNPWMVPS LNPPPVNKGK ANIFYKTRMC AKFRAGTCRN GELCNFAHGI EDLRQPPSNW QEIVGPPAG QDRERERERE RERERPSLAP VVNNNWEDDQ KIILRMKLCR KFCFGEECPY GDRCNFIHED LSKFREDGSK LRESSVISVG ATAADQPSDT ASNLIEVNRQ GSIPVPAPMN NGGVVKTYYW KTRLCMKFDI TGQCPFGDKC HFAHGQAEHL NSVGRVEGEA MNAVASVKNQ AVVPANEAFA MKPITQVTAD SSGLNEEGRR KKCLLKWSDS KKINRIYGDW IDDLPGVQKS TKPVES
Source	Mammalian cell
Target Names	At3g19360
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 39 Short name= AtC3H39
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.