



Recombinant Arabidopsis thaliana High mobility group B protein 10 (HMGB10)

Product Code	CSB-EP888817DOA-B
Abbreviation	HMGB10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LTT3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSTDISPPYS QTHVEPVNGY PSDNKRCDDS SVPAKYDDL V RNSALFWEKL RAFLGLTSKT LKVPTVGGNT LDHRLFIEV TSRGGIERVV KDRKWKEVIG AFSFP TTITS ASFVLRKY YL KFLFQLEHVY YLEKPVSS LQ STDEALKSLA NESPNPEEGI DEPQVGYEVQ GFIDGKFDSG YLVTMKLGSQ ELKGVLYHIP QTPSQSQQTM ETPSAIVQSS QRRHRKSKL AVVDTQKPKC HRSGYNFFFA EQYARLKPEY HGQERSITKK IGHMWSNLTE SEKQVYQDKG VKDVERYRIE MLEYKSSHES GATASTVAQ
Source	E.coli
Target Names	HMGB10
Protein Names	Recommended name: High mobility group B protein 10 Alternative name(s): Nucleosome/chromatin assembly factor group D 10
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.