



Recombinant Arabidopsis thaliana F-box protein At3g17710 (At3g17710)

Product Code	CSB-EP888827DOA
Abbreviation	At3g17710
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LUM6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASVKLPWDL EEEILSRLPP RSLVRFRTVC KHWNGLFSDK RFVKKHLVRA RPQFIFLTES KKMYSIEIDL GGTIEVREVP YDFHCQPMKK NFFTIMACDG LLFRDFWKQG VAVWNPWLRQ VGWIEYEDKG FRFCGVGYDS CKPDKCYKIL GYFNCTRRLS DSLQEGYYQA AIYECASQAF KFIDTPNPFN LWPNKDPLSV NGNLYWLAHN HPETLEYFIE TFDIFSMEIFK PFCLLPCRKD FGSNELVLAV FKEDRFSLK QCFETTKIEI WVTMKMIDRE EEVWIKFMT LPTTNLPNLD DTYCCSSYFI FDKTIIMCCG DYKTGATCIY MVRGDMFKKI QINSGIFQFS YCVYLPNLIS VPFESHQV
Source	E.coli
Target Names	At3g17710
Protein Names	Recommended name: F-box protein At3g17710
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.