



Recombinant Arabidopsis thaliana Putative CCR4-associated factor 1 homolog 8 (CAF1-8)

Product Code	CSB-MP888855DOA
Abbreviation	CAF1-8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LXM4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSLIEDCLRS YRFIAIDTEF PSTLRETTQH ATDEERYMDM SFSVDRAKLI QLGLTLFDIN GRIGGTWEIN FSDFGVDDAR NEKSIEFLRR NGLDLRKIRE EGIRIEGFFS EMFWMLKKTR RNITWVTFHG SYDIAYLLKG FTGEALPVT ERFSKAVARV LGSVYDLKVM AGRCEGLSSR LGLETLAHEF GLNRVGTAAH AGSNNELTAM VFAKVLSPFP LFLRFGSLEL RTIESEAIV
Source	Mammalian cell
Target Names	CAF1-8
Protein Names	Recommended name: Putative CCR4-associated factor 1 homolog 8 EC=3.1.13.4
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.