



# Recombinant Arabidopsis thaliana 6-phosphofructokinase 6 (PFK6)

<b>Product Code</b>	CSB-EP888881DOA-B
<b>Abbreviation</b>	PFK6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9M076
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASNGVDEQI KLVEGPAGYV LEDVPHLSY ILDLPTYPNP LQSNAAYSVV RQYFVDEDDT VQEKIVVHKD SPRGTHFRRA GPRQKVYFKP SDVRACIVTC GGLCPGLNTV IREIVCGLHF MYGVTEVIGV DCGFRGFYSK NTVALTPKTV SDIHKRGGTI LGTSRGGHDT SKIVDNIQDR EINQVYIIGG DGTQKGANAI YKEIRRRGLK VAVAGIPKTI DNDIPVIDKS FGFDTAVEEA QRAINAAHVE ATSVENGIGI VKLMGRYSGF IAMYATLASR DVDCCLIPES PFYLEGKGGI YEFIKRLRE NGHMVIVIAE GAGQDLVAES IEQQDASGNK LLKDVGLWMS LKIKEYFAKH NVMDITLKYI DPTYMIRAIP ANASDNVYST LLAQSAVHGA MAGYTGFBVSG LVNGRHTYIP FNRITERQNK VVITDRMWAR MLSSTNQPSF MNPPKGTTEF TD
<b>Source</b>	E.coli
<b>Target Names</b>	PFK6
<b>Protein Names</b>	Recommended name: 6-phosphofructokinase 6 EC= 2.7.1.11 Alternative name(s): Phosphofructokinase 6 Phosphohexokinase 6
<b>Expression Region</b>	1-462
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.