



Recombinant Arabidopsis thaliana 6-phosphofructokinase 6 (PFK6)

Product Code	CSB-MP888881DOA
Abbreviation	PFK6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M076
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASNGVDEQI KLVEGPAGYV LEDVPHLSDY ILDLPTYPNP LQSNAAYSVV RQYFVDEDDT VQEKIVVHKD SPRGTHFRRA GPRQKVYFKP SDVRACIVTC GGLCPGLNTV IREIVCGLHF MYGVTEVIGV DCGFRGFYSK NTVALTPKTV SDIHKRGGTI LGTSRGGHDT SKIVDNIQDR EINQVYIIGG DGTQKGANAI YKEIRRRGLK VAVAGIPKTI DNDIPVIDKS FGFDTAVEEA QRAINAAHVE ATSVENGIGI VKLMGRYSGF IAMYATLASR DVDCCLIPES PFYLEGKGGI YEFIAKRLRE NGHMVIVIAE GAGQDLVAES IEQQDASGNK LLKDVGLWMS LKIKEYFAKH NVMDITLKYI DPTYMIRAIP ANASDNVYST LLAQSAVHGA MAGYTGFBVSG LVNGRHTYIP FNRITERQNK VVITDRMWAR MLSSTNQPSF MNPPKGTTEF TD
Source	Mammalian cell
Target Names	PFK6
Protein Names	Recommended name: 6-phosphofructokinase 6 EC= 2.7.1.11 Alternative name(s): Phosphofructokinase 6 Phosphohexokinase 6
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.