



Recombinant Arabidopsis thaliana F-box/kelch-repeat protein SKIP20 (SKIP20)

Product Code	CSB-YP888897DOA
Abbreviation	SKIP20
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M1Y1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGVSKKKSGE IRGDLIPGLP EELAIECLVR VPFQFHSSIK SVCRSWKCVI SSRSFIKERI GFGKAESLLC LVQPLTSPPS PAMMEGGEMS QKKKEEEEGE SQMTQQLLQP RITGTPLYGL SVYNATLDTW HRVAIPERIP LFCECVAIQD AGKVLLIGGW DPETLQPVRD VFVLDFPAGE GSGRRFRRGR PMSAARSFFA CASVGSTKVY VAGGHDDQKN ALRSAEVYDV EKDEWSMLPP MTEGRDECHG FSMATDPGFC VLSGYGTETQ GQFRSDGEIY DPITNSWSTI ENVWFPDLS PRGRTAAAAA EFGDFRGCR LWCDFIDSERQ SQPHWEVEDD SMKWKVIMDT IRLPVTTMTS VFAGSLSGQA VAMIGGGGEE SGTMMVKTTA EKNGGKWSHV NTPSGFSSLP FSCSSIYV
Source	Yeast
Target Names	SKIP20
Protein Names	Recommended name: F-box/kelch-repeat protein SKIP20 Alternative name(s): SKP1-interacting partner 20
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.