



Recombinant Arabidopsis thaliana Probable plastid-lipid-associated protein 9, chloroplastic (PAP9)

Product Code	CSB-EP888903DOA
Abbreviation	PAP9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M2P7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	AMVQDSV QGIPSVYARE MERLSAKESL ILAFNDAGGF EALVTGKITD MQKIDVNERI TNLERLNPTP RPTTSPYLEG RWSFEWFGVN TPGSLAVRVM FERFPSTLVS LSNMEIFIKD NNTKATANIK LLNSIENKIT LSSKLTIEGP LRMKEEYLEG LLESPTVIEE AVPDQLRGLL GQATTTLQQL PEPIKDTLAN GLRIPLGGTY QRFFMISYLD DEILIVRDTA GVPEVLTRVE TSSPMSSSSV VENLEYNS
Source	E.coli
Target Names	PAP9
Protein Names	Recommended name: Probable plastid-lipid-associated protein 9, chloroplastic Alternative name(s): Fibrillin-9 Short name= AtPGL34
Expression Region	54-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.