



# Recombinant Arabidopsis thaliana B3 domain-containing protein REM21 (REM21)

<b>Product Code</b>	CSB-EP888931DOA-B
<b>Abbreviation</b>	REM21
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9M8K2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MWATKSRRCA VHTSSMASGD VLPRFFTVFL SHCSSESMVI PRSYYNLLPR PLPKTAILIG TGGRFWKVAM TSKQEQQVYFE QGWGNFVADN QLKEGEFLTF VFDGHKSYEV SIYGRGDCKE TRAVIQVEEI SDDTEDDNVS LHSPSNVSLD SLSNDSHHST SNVSLRSLSN DSLHGDAEIE SDSEYSPENL PSASISVESV EVVNPTTSRQ RSYKRKTIEN PHLYLDDPNN VCFETCLKLR KFELLVHAQL VKDYGLIFSD NVDYIDGYGK LTAKTTKWAD QRVCINKWQK ICERNQFTEN DSILCEILRN EDKVVYAIKI HIFRDAAAST
<b>Source</b>	E.coli
<b>Target Names</b>	REM21
<b>Protein Names</b>	Recommended name: B3 domain-containing protein REM21 Alternative name(s): Protein REPRODUCTIVE MERISTEM 21
<b>Expression Region</b>	1-330
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.