



# Recombinant Arabidopsis thaliana F-box/LRR-repeat protein At3g03030 (At3g03030)

<b>Product Code</b>	CSB-EP888933DOA
<b>Abbreviation</b>	At3g03030
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9M8U4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDLVSSLPDD VRCHILSFLT TKESALTSVL SKKWRNLFTL VPNLDDLDSE YLHPEETKWE REEILDSFLD FVERVLSLQG NGPIRRFSLR CESGVPPARL NRWLCKVLQR GVADLELIID FEDGYLLPRE LRVSETLVNL RLKSEFGCIH WWPGAEGTFL PRLQSLDISS AMVFCDDKLQ KLLPCFPVIE ELRLVGMIEWI DSHVTVSCAT LTNQLLLSST GSRTSKNPKS VSFDTPNLLS LAYS DLVAED YPLVNMKSLF KARIILGVDD DQIARMRAPT NYLSEDDVVR RFGNVVKLMK GIQNVQQLHL HSDTLEVLSM CCESMPVFNN VKTLGIKSDY DRGWQAVPAL LRNCPHLETL TFEGLVHHVT DKCGDACDCI FRKDKGRSLK SCPVKVVEIK EFRVTMKEMH LIEHFLDNFP CLKEMKIYRE EYDPSFPGND PQVDDLILDM IEEYKDIYSC NVQLLEPCDK WDEQ
<b>Source</b>	E.coli
<b>Target Names</b>	At3g03030
<b>Protein Names</b>	Recommended name: F-box/LRR-repeat protein At3g03030
<b>Expression Region</b>	1-474
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.