



Recombinant Arabidopsis thaliana F-box/LRR-repeat protein At3g03030 (At3g03030)

Product Code	CSB-MP888933DOA
Abbreviation	At3g03030
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M8U4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDLVSSLPDD VRCHILSFLT TKESALTSVL SKKWRNLFTL VPNLDDLDSE YLHPEETKWE REEILDSFLD FVERVLSLQG NGPIRRFSLR CESGVPPARL NRWLCKVLQR GVADLELIID FEDGYLLPRE LRVSETLVNL RLKSEFGCIH WWPGAEGTFL PRLQSLDISS AMVFCDDKLQ KLLPCFPVIE ELRLVGMIEWI DSHVTVSCAT LTNQLLLSST GSRTSKNPKS VSFDTPNLLS LAYS DLVAED YPLVNMKSLF KARIILGVDD DQIARMRAPT NYLSEDDVVR RFGNVVKLMK GIQNVQQLHL HSDTLEVLSM CCESMPVFNN VKTLGIKSDY DRGWQAVPAL LRNCPHLETL TFEGLVHHVT DKCGDACDCI FRKDKGRSLK SCPVKVVEIK EFRVTMKEMH LIEHFLDNFP CLKEMKIYRE EYDPSPFGND PQVDDLILDM IEEYKDIYSC NVQLLEPCDK WDEQ
Source	Mammalian cell
Target Names	At3g03030
Protein Names	Recommended name: F-box/LRR-repeat protein At3g03030
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.