



Recombinant Arabidopsis thaliana Ethylene-responsive transcription factor ERF117 (ERF117)

Product Code	CSB-BP888935DOA
Abbreviation	ERF117
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M9B2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MISFREENID LNLIKTISVI CNDPDATDSS SDDESISGNN PRRQIKPKPP KRYVSKICVP TLIKRYENVV NSTGNKAAGN RKTSSGFKGV RRRPWGKFAA EIRNPFKKR KWLGTPTTEE EAAEAYQKSK REFDERLGLV KQEKDLVDLT KPCGVRKPEE KEVTEKSNCK KVNKRIVTDQ KPFGCGYNAD HEEEGVISKM LEDPLMTSSI ADIFGDSAVE ANDIWVDYNS VEFISIVDDF KFDFVENDRV GKEKTFGFKI GDHTKVNQHA KIVSTNGDLF VDDLLDFDPL IDDFKLEDFP MDDLGLLGDP EDDDFSWFNG TTDWIDKFL
Source	Baculovirus
Target Names	ERF117
Protein Names	Recommended name: Ethylene-responsive transcription factor ERF117
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.