



Recombinant Arabidopsis thaliana Interactor of constitutive active ROPs 4 (ICR4)

Product Code	CSB-EP888937DOA-B
Abbreviation	ICR4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M9F9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MPKPSIRGSE LPQRQSPRLR TLLSTSSDP HHLSRPITDR SPKLGLDRRS PRSGGPHTDP LSQKKLGSRI SGLESQLGQA QEELRLLKQQ LAKAEAAKKR AQEELHRKKS KKPNTPAPER DDIPGDGHQE TDVFEVLDEK AKESEKTKND ELASKEDQIN VLKARLYDLE KERVSLSEEN ETLKDQLKKT DTEMSCAKAK EDEIASKVSQ IGEELEESNE TTAKLKKKLE SVEEAKETLE AEMKKLKVQT EQWRKAADAA AAVLSGGVEM NGRFSEQCGS MEKHFAGRFV GSPGMADDSD DGSGKRKSSG KKMFGDLWKK KGQK
Source	E.coli
Target Names	ICR4
Protein Names	Recommended name: Interactor of constitutive active ROPs 4 Alternative name(s): ROP interactive partner 2
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.