



Recombinant Protein SEC13 homolog (npp-20)

Product Code	CSB-EP889030CXY-B
Abbreviation	npp-20
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9N4A7
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTTVRQRIDT QHRDAIHDAQ LNIYGSRLAT CGSDRLVKIF EVRPNGQSYP MAELVGHSGP VWKVSWAHPK YGGLLASASY DKKVIWNEQ QGRWQKAYEW AAHEASTTCV AFAPHQYGLM LASASADGDI GILRYDNSSN EWISSKIQKC HEQGVNSVCW APGSADPAAK KRLVSAGNDK NVKIWAFDDA TNEWILEKTL AGHTDFVREA AWCPVTNNGQ HTIVSCGMEG NLVLFRTSNI ETEEWKAKLL ETAPCALYHS SFSPCGSFLS VAGDDNVITI WRENLQGQWI KVPRDNKERE GMSQAVGAPG AQR
Source	E.coli
Target Names	npp-20
Protein Names	Recommended name: Protein SEC13 homolog Short name= CeSEH13R Alternative name(s): Nuclear pore complex protein 20
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.