



Recombinant Putative T-box protein 39 (tbx-39)

Product Code	CSB-MP889036CXY
Abbreviation	tbx-39
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9NA56
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	<p>MFSQAQFNVT MAEEDRWKQW YPNMEMIVNN KGPKLIFPEL NFNITGLEEK SRYVVLLSIE KYDNIRYGYR NGKWGPSKVR HATKKEQEIK YFLHPDGTKL GKELMKETIK FDTVIRITNHK KFMDKDNVFF VETMHKYVPV LTVKNITNVE SANHSMRMEV AQFFPVTVYN QESIGNWWSK FHKNSTYGNR LDGGNKRKNT DSSEERTSKR SKNETEIAVS DILQVASQSE NYNESTNSGR LQNEISSIH QFPSTSYQSQ YPHAYPTVNT PPIYVQQFTA LFDEIHNQFD PKTKQNRVAK NVQYLSSVPQ FAINQKENIH QNAPDYPVSQ FPNQYPYGSY PYYYPYFQSS QQFPHLNYSM NNSIYSDSSF QSSFSAENSY FDENGNPIHY PYYSSN</p>
Source	Mammalian cell
Target Names	tbx-39
Protein Names	Recommended name: Putative T-box protein 39
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>