



Recombinant Cell death abnormality protein 2 (ced-2)

Product Code	CSB-YP889052CXY
Abbreviation	ced-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NHC3
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTTNGFDPFE WRSFYFPGMS REEAHKLLGE PQVSIQTFLM RDSSRPGEYS LTVREADEGN AVCHYLIERG EPKEDGTAAA GVKIANQSFP DIPALLNHFK MRVLTEASLL AAYKKPIIEV VVGTFKFTGE RETDLPFEQG ERLEILSKTN QDWWEARNAL GTTGLVPANY VQIQMEFHND RTSKGASQSS IGSSGGGAER FSSASTSSDN IELQPRLPAK AKVTFDRVPN AYDPTQLRVK KGQTVLVTQK MSNGMYKAEL DGQIGSVPHY YLRFTAVSE
Source	Yeast
Target Names	ced-2
Protein Names	Recommended name: Cell death abnormality protein 2 Alternative name(s): Cell-corpse engulfment protein CED-2
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.