



Recombinant Human Exosome complex component RRP41 (EXOSC4)

Product Code	CSB-MP889069HU
Abbreviation	EXOSC4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NPD3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AGLELLSDQ GYRVDGRRAG ELRKIQARMG VFAQADGSAY IEQGNTKALA VYGPHEIRG SRARALPDRA LVNCQYSSAT FSTGERKRRP HGDRKSCEMG LQLRQTFEAA ILTQLHPRSQ IDIYVQVLQA DGGTYAACVN AATLAVLDAG IPMRDFVCAC SAGFVDGTAL ADLSHVEEAA GGPQLALALL PASGQIALLE MDARLHEDHL ERVLEAAAQA ARDVHTLLDR VVRQHVREAS ILLGD
Source	Mammalian cell
Target Names	EXOSC4
Protein Names	Recommended name: Exosome complex component RRP41 Alternative name(s): Exosome component 4 Ribosomal RNA-processing protein 41 p12A
Expression Region	2-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.