



Recombinant Human Lysophosphatidic acid phosphatase type 6 (ACP6)

Product Code	CSB-EP889072HU
Abbreviation	ACP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NPH0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ELQEADGQ CPVDRSLLKL KMQVQVFRHG ARSPLKPLPL EEQVEWNPQL LEVPPQTQFD YTVTNLAGGP KPYSYDSQY HETTLKGGMF AGQLTKVGMQ QMFALGERLR KNYVEDIPFL SPTFNPQEVF IRSTNIFRNL ESTRCLLAGL FQCQKEGP II IHTDEADSEV LYPNYQSCWS LRQRTRGRRQ TASLQPGISE DLKKVKDRMG IDSSDKVDFE ILLDNVAAEQ AHNLPSCPML KRFARMIEQR AVDTSLYILP KEDRESLQMA VGPFLHILES NLLKAVDSAT APDKIRKLYL YAAHDVTFIP LLMTLGIFDH KWPPFAVDLT MELYQHLESK EWFVQLYYHG KEQVPRGCPD GLCPLDMFLN AMSVYTLSPE KYHALCSQTQ VMEVGNEE
Source	E.coli
Target Names	ACP6
Protein Names	Recommended name: Lysophosphatidic acid phosphatase type 6 EC= 3.1.3.2 Alternative name(s): Acid phosphatase 6, lysophosphatidic Acid phosphatase-like protein 1 PACPL1
Expression Region	33-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.