



Recombinant Human PR domain zinc finger protein 8 (PRDM8)

Product Code	CSB-YP889084HU
Abbreviation	PRDM8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NQV8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEDTGIQRGI WDGDAKAVQQ CLTDIFTSVY TTCDIPENAI FGPCVLSHTS LYDSIAFIAL KSTDKRTVPY IFRVDTSAAN GSSEGLMWLR LVQSARDKEE QNLEAYIKNG QLFYRSLRRI AKDEELLVWY GKELTELLLL CPSRSHNKMN GSSPYTCLEC SQRQFQFEFPY VAHLRFRCPK RLHSADISPO DEQGGGVGK DHGGGGGGGK DQQQQQQEAP LGPGPKFCKA GPLHHYSPS PESSNPSAAA GGSSAKPSTD FHNLARELEN SRGGSSCSPA QSLSSGSGSG GGGGHQEDEL SPDGIATGGG KGKRKFPEEA AEGGGGAGLV GGRGRFVERP LPASKEDLVC TPQQYRASGS YFGLEENGRL FAPPSPETGE AKRSFVEVK KAARAASLQE EGTADGAGVA SEDQDAGGGG GSSTPAAASP VGAEKLLAPR PGGPLPSRLE GGSPARGSAF TSVPQLGSAG STSGGGGTGA GAAGGAGGGQ GAASDERKSA FSQPARSFSQ LSPLVLGQKL GALEPCHPAD GVGPTRLYP AADPLAVKLQ GAADLNGGCG SLPSGGGGLP KQSPFLYATA FWPKSSAAAA AAAAAAAGP LQLQLPSALT LLPPSFTSLC LPAQNWCAKC NASFRMTSDL VYHMRS HHK EYAMEPLVKR RREEKLCPI CNESFRERHH LSRHMTSHN
Source	Yeast
Target Names	PRDM8
Protein Names	Recommended name: PR domain zinc finger protein 8 Alternative name(s): PR domain-containing protein 8
Expression Region	1-689
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.