



# Recombinant Human Something about silencing protein 10 (UTP3)

<b>Product Code</b>	CSB-EP889085HU-B
<b>Abbreviation</b>	UTP3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NQZ2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVGRSRRRGA AKWAAVRAKA GPTLTDENGD DLGLPPSPGD TSYYQDQVDD FHEARSRAAL AKGWNEVQSG DEEDGEEEE EVLALDMDDE DDEDGGNAGE EEEEEENADDD GGSSVQSEAE ASVDPSLSWG QRKKLYYDTD YGSKSRGRQS QQEAEEEEERE EEEEEAQIIQR RLAQUALQEDD FGVAVVEAFA KPVPQVDEAE TRVVKDLAKV SVKEKLMMLR KESPELLELI EDLKVKLTEV KDELEPILLEL VEQGIIPPGK GSQYLRTKYN LYLNYCSNIS FYLILKARRV PAHGHPVIER LVTYRNLINK LSVVDQKLSS EIRHLLTLKD DAVKKELIPK AKSTKPKPKS VSKTSAAACA VTDLSDDSD DEKAKLKYYK EIEDRQKLKR KKEENSTEEQ ALEDQNAKRA ITYQIAKNRG LTPRRKKIDR NPRVKHREKF RRAKIRRRGQ VREVRKEEQR YSGELSGIRA GVKKSILK
<b>Source</b>	E.coli
<b>Target Names</b>	UTP3
<b>Protein Names</b>	Recommended name: Something about silencing protein 10 Alternative name(s): Charged amino acid-rich leucine zipper 1 Short name= CRL1 Disrupter of silencing SAS10 UTP3 homolog
<b>Expression Region</b>	1-479
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.