



Recombinant Human Something about silencing protein 10 (UTP3)

Product Code	CSB-EP889085HU
Abbreviation	UTP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NQZ2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVGRSRRRGA AKWAAVRAKA GPTLTDENGD DLGLPPSPGD TSYYQDQVDD FHEARSRAAL AKGWNEVQSG DEEDGEEEE EVLALDMDDE DDEDGGNAGE EEEEEENADDD GGSSVQSEAE ASVDPSLSWG QRKKLYYDTD YGSKSRGRQS QQEAEEEEERE EEEEEAQIIQR RLAQALQEDD FGVAVVEAFA KPVPQVDEAE TRVVKDLAKV SVKEKLMMLR KESPELLELI EDLKVKLTEV KDELEP LLEL VEQGIIPPGK GSQYLRTKYN LYLNYCSNIS FYLILKARRV PAHGHPVIER LVTYRNLINK LSVVDQKLSS EIRHLLTLKD DAVKKELIPK AKSTKPKPKS VSKTSAACA VTDLSDDSD DEKAKLKYYK EIEDRQKLKR KKEENSTEEQ ALEDQNAKRA ITYQIAKNRG LTPRRKKIDR NPRVKHREKF RRAKIRRRGQ VREVRKEEQR YSGELSGIRA GVKKSILK
Source	E.coli
Target Names	UTP3
Protein Names	Recommended name: Something about silencing protein 10 Alternative name(s): Charged amino acid-rich leucine zipper 1 Short name= CRL1 Disrupter of silencing SAS10 UTP3 homolog
Expression Region	1-479
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.