



Recombinant Human Phosphoribosyltransferase domain-containing protein 1 (PRTFDC1)

Product Code	CSB-MP889092HU
Abbreviation	PRTFDC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NRG1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AGSSEEAPD YGRGVVIMDD WPGYDLNLFT YPQHYYGDLE YVLIPHGIIV DRIERLAKDI MKDIGYSDIM VLCVLKGGYK FCADLVEHLK NISRNSDRFV SMKVDFIRLK SYRNDQSMGE MQIIGDDLS TLAGKNVLIV EDVVGTRTM KALLSNIKEY KPNMIKVASL LVKRTSRSDG FRPDYAGFEI PNLFVVGAL DYNEYFRDLN HICVINEHGK EKYRV
Source	Mammalian cell
Target Names	PRTFDC1
Protein Names	Recommended name: Phosphoribosyltransferase domain-containing protein 1
Expression Region	2-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.