



Recombinant Human E3 ubiquitin-protein ligase RNF146 (RNF146)

Product Code	CSB-YP889110HU
Abbreviation	RNF146
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NTX7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMAGCGEIDH SINMLPTNRK ANESCSNTAP SLTVPECAIC LQTCVHPVSL PCKHVFCYLC VKGASWLGKR CALCRQEIPE DFLDKPTLLS PEELKAASRG NGEYAWYYEG RNGWWQYDER TSRELEDAFS KGKKNTEMLI AGFLYVADLE NMVQYRRNEH GRRRKIKRDI IDIPKKGVAG LRLDCDANTV NLARESSADG ADSVSAQSGA SVQPLVSSVR PLTSVDGQLT SPATPSPDAS TSLEDSFAHL QLSGDNTAER SHRGEGEEDH ESPSSGRVPA PDTSIEETES DASSDSEDVS AVVAQHSLTQ QRLLVSNANQ TVPDRSDRSG TDRSVAGGGT VSVSVRSRRP DGQCTVTEV
Source	Yeast
Target Names	RNF146
Protein Names	Recommended name: E3 ubiquitin-protein ligase RNF146 EC= 6.3.2.- Alternative name(s): Dactylidin RING finger protein 146
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.