



# Recombinant Human ATP-binding cassette sub-family F member 3 (ABCF3)

<b>Product Code</b>	CSB-EP889115HU
<b>Abbreviation</b>	ABCF3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NUQ8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ATCAEILRS EFPEIDGQVF DYVTGVLHSG SADFESVDDL VEAVGELLQE VSGDSKDDAG IRAVCQRMYN TLRLAEPQSQ GNSQVLLDAP IQLSKITENY DCGTKLPGLL KREQSSTVNA KKLEKAEARL KAKQEKRSEK DTLKTSNPLV LEEASASQAG SRKESRLESS GKNKSYDVRI ENFDVSFGDR VLLAGADVNL AWGRRYGLVG RNGLGKTTLL KMLATRSLRV PAHISLLHVE QEVAGDDTPA LQSVLESDSV REDLLRRERE LTAQIAAGRA EGSEAAELAE IYAKLEEIEA DKAPARASVI LAGLGFTPCKM QQQPTREFSG GWRMRLALAR ALFARPDLLL LDEPTNMLDV RAILWLENYL QTWPSTILVV SHDRNFLNAI ATDIIHLHSQ RLDGYRGDFE TFIKSKQERL LNQQREYEAQ QQYRQHIQVF IDRFRYNANR ASQVQSKLKM LEKLPKPV DKESEVVMKF PDGFEKFSPP ILQLDEVDFY YDPKHVIFSR LSVSADLESR ICVVGENGAG KSTMLKLLL DLAPVRGIRH AHRNLKIGYF SQHHVEQLDL NVSAVELLAR KFPGRPEEEY RHQLGRYGIS GELAMRPLAS LSGGQKSRVA FAQMTMPCPN FYILDEPTNH LDMETIEALG RALNNFRGGV ILVSHDERFI RLVCRELWVC EGGGVTRVEG GFDQYRALLQ EQFRREGFL
<b>Source</b>	E.coli
<b>Target Names</b>	ABCF3
<b>Protein Names</b>	Recommended name: ATP-binding cassette sub-family F member 3
<b>Expression Region</b>	2-709
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.